



# NOBLE MANES CLINIC

SERVICES FOR BREEDERS BY BREEDERS  
Animal Clinic, Diagnostic Laboratory Testing, Breeder Support & Supplies



## Toxicology Screen

### IDEXX TEST CODE

IDEXX 1433

### IDEXX TEST TYPE

Gas Chromatography Mass Spectrometry (GC-MS)

### REQUIRED SUPPLIES

Collection Pack: Universal

### TURNAROUND TIME

15 - 23 days

### TESTS INCLUDED

This is an extensive screen for a wide variety of pesticides, drugs, pollutants, and industrial chemicals.

### SPECIMEN REQUIREMENTS: 1433

10 mL urine divided into two sterile container, WTT (White Top Tube) preferred; keep refrigerated. Please specify collection method. Optionally: 100 g rumen/crop/stomach contents in a sterile container, RTFC (Red Top Fecal Container) preferred; keep refrigerated. This test can also be run on the suspected material: 400 g in a sterile container, RTFC (Red Top Fecal Container) preferred; keep refrigerated.

### STORAGE AND STABILITY

Urine up to 24 hours at 36° to 46° F (2° to 8° C), keep refrigerated. Stomach Contents up to 24 hours at 36° to 46° F (2° to 8° C), keep refrigerated.

### TEST COMPONENTS: REFERRED

- Toxicology Urine Screening

### NOTES

Submission of the most appropriate sample type, based on clinical history and toxins suspected, is important for best results. If there is uncertainty about which sample type to submit, consultation with a clinical toxicologist is recommended.

### INTERFERENCE

[Mass-Spec] poor sample quality / aged sample / exposure to temperature extremes / formalin / saline.