



NOBLE MANES CLINIC

SERVICES FOR BREEDERS BY BREEDERS

Animal Clinic, Diagnostic Laboratory Testing, Breeder Support & Supplies



Urinalysis

IDEXX TEST CODE

IDEXX 910

IDEXX TEST TYPE

Chemistry

REQUIRED SUPPLIES

Collection Pack: Universal

TURNAROUND TIME

8:00am for Urinalysis

TESTS INCLUDED

Urinalysis: physical, chemical and microscopic analyses

SPECIMEN REQUIREMENTS: 910

5 mL urine in a sterile container, WTT (White Top Tube) preferred; keep refrigerated. Please specify collection method.

STORAGE AND STABILITY

Urine up to 24 hours at 36° to 46° F (2° to 8° C), keep refrigerated.

TEST COMPONENTS: REFERRED

- Bacteria
- Bilirubin
- Blood
- Casts
- Clarity
- Color
- Crystals
- Cystatin B
- Epi Cell
- Glucose
- Ketones
- Mucus
- Ph
- Protein
- RBC
- Specific Gravity

- Urobilinogen
- WBC

NOTES

Do not submit syringes with needles.

INTERFERENCE

[Urinalysis] poor sample quality / aged sample / clot activator / EDTA / myoglobinuria or hemoglobinuria may cause a positive for blood on dipstick in the absence of red blood cells on sediment / captopril and ace-inhibitors may cause false negative for blood on dipstick.